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Certificado de calidad
Quality certification
Bionova® BT10

Esterilización por OE / EO sterilization
Bacillus atrophaeus ATCC 9372

Población / _____ UFC /
Population _____ CFU

Valor D / _____ min.

(54 °C, 60% HR, OE 600 mg/L)

D - value

(54 °C, 60% RH, EO 600 mg/L)

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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Parameters determined at time of manufacture according to ISO 11138-1: 2006, ISO 11138-2:2006 and IRAM 37102: 1999 (Parts 1 and 2) standards. The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant.

ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Composition

Each tube contains a population of spores of *Bacillus atrophaeus* ATCC 9372 soaked on a carrier. It also has a growth indicator medium contained within the glass ampoule.

Product description

Bionova® BT10 Biological Indicator is specifically designed for monitoring Ethylene Oxide (EO) sterilization processes. If the sterilization process has not been successful, the indicator medium will turn to yellow after incubation at 37 ± 2 °C, thus indicating the presence of living *Bacillus atrophaeus* spores.

If the sterilization process has been successful, the indicator medium will remain the original color after incubation. The final readout should be performed after 48 hours of incubation at 37 ± 2 °C.

Warning: Do not use Bionova® BT10 Biological Indicator for monitoring Steam sterilization cycles, Dry Heat, Chemical Vapor, Radiation or other sterilization processes. Do not re use biological indicators.

Storage

Best stored in the original box under the following conditions: temperature between 10-30 °C, 30-80 % Relative Humidity.

Do not freeze.

Do not store biological indicators near sterilizing agents or other chemical products.

Instructions

1. Identify the Bionova® BT10 Indicator tube by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.
2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for the sterilizing agent (EO). Generally a problematic area is the center of the load.

3. Sterilize as usual.

4. After the sterilization process has finished, you can:

A. Remove the biological indicator from the sterilization package for incubation prior to material aeration.

B. Aerate the package containing the biological indicator with the rest of the load and then remove the biological indicator for incubation. **CAUTION!** Wear safety glasses and gloves when removing the Bionova® BT10 Biological Indicator from the sterilized package.

WARNING! Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5. Check the chemical indicator on the label of the biological indicator. A color change to green confirms that the biological indicator has been exposed to Ethylene Oxide. **IMPORTANT:** This color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

6. Crush the ampoule contained in the biological indicator.

This may be done in 3 different ways:

A- With an ampoule crusher.

B- Manually. To do this, take the tube by putting your index finger and thumb on the space between the cap and the line of chemical indicator and then compress the tube.

C- With the ampoule crusher placed within the top of the incubator's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator.

IMPORTANT: Use a non-sterilized biological indicator as a positive control in order to ensure that correct incubation conditions were met, capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals. Both the positive control indicator and the processed indicator should belong to the same batch.

7. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 48 hours at 37 ± 2 °C. Readout should be performed at convenient intervals of 10 hours. A color change to yellow of the growth indicator medium means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained (the sterilization process was satisfactory). The positive control indicator should show a yellow color change for the results to be valid. Record the positive results and discard them immediately as shown below.

WARNING! Do not reuse the sterilizer until the biological indicator test results are negative (processed indicator remains the original color).

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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